# kliniderm® alginate standard & alginate cavity

## INSTRUCTIONS FOR USE

### ENGLISH

**Product Description**
Kliniderm® Alginate dressings are sterile, nonwoven calcium/sodium alginates consisting of a high G (Guluronic acid) alginate. Kliniderm® Alginate dressings are soft, clear gelling, highly absorbent and ensure intact removal. It used correctly does not stick to the wound.

Not made with natural rubber latex. Medical device.

**Indications**
Kliniderm® Alginate dressings are indicated for moderate to heavily exuding chronic and acute wounds including pressure ulcers, leg ulcers, diabetic ulcers, post-operative wounds and donor sites. Kliniderm® Alginate dressings can also be used to control minor bleeding in superficial wounds.

**Contra-indications**
Kliniderm® Alginate dressings are not indicated to control heavy bleeding or for use as surgical sponge.

### Directions for use

Application:


- Select a size of Kliniderm® Alginate dressing that is slightly larger than the wound. The dressing can be cut to size if necessary.
- For cavity wounds, loosely fill ensuring that the dressing does not overlap the wound margins.
- Cover and secure with a semi occlusive secondary dressing, for example Kliniderm® Film or Kliniderm® Foam dressing.

**Dressing Change and Removal:**

- Kliniderm® Alginate dressings can remain in place for up to 7 days, dependent on patient condition and the level of exudate.
- Replace Kliniderm® Alginate when the dressing has reached its absorbent capacity or whenever good wound care practice dictates that the dressing should be changed.
- Gently remove the dressing from the wound bed and dispose according to local procedures and guidelines. If removal appears to be difficult, saturate the dressing with sterile saline solution.
- If required, irrigate the wound site according to good wound care practice.

### Precautions

- In case of signs of infection, consult a healthcare professional for adequate infection treatment.
- Do not use on patients with a known sensitivity to alginates.
- Sterile if pouch and seal are intact; do not use if damaged.
- Do not re-use. Re-use can cause serious infections.
- Do not resterilize.
- For external use only.
- For use under the guidance of a healthcare professional.

In case a serious incident has occurred in relation to the dressing inform Medico BV as listed on the packaging, your local distributor and your local healthcare authority.

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