

Devon™ Positioning Products

Information for patients

What are Devon™ Positioning Products?

The range includes Devon foot and heel and Devon utility pad. They are designed to reduce the risk of pressure ulcers and nerve damage while maintaining proper circulation, to be used as part of a pressure ulcer prevention or treatment strategy.

Where can they be used?

They can be used on any area of the body that requires support stability and protection for the legs, ankles, heels, feet as well as elbows, knees and head.

Which way should they be applied?


The bumpy (convoluted) side should be placed against the skin.

How to apply?

Your Healthcare Professional will advise on how to apply these products.


Devon foot and heel

1. Position foot inside ensuring that the seam is not below the resting area of the foot
2. Fix in place using the Velcro® strap. Fit needs to be comfortable without being too tight



Devon utility pad

1. Position underneath the area of the body that requires support, stability and protection – bumpy side up



How often should these products be changed?

These can be used until they are no longer needed. Your Healthcare Professional should replace these products if they become damaged or dirty.

Can a dressing still be worn?

Yes, if you have a wound it is important that it is covered with an appropriate dressing, your Healthcare Professional will provide advice on this.

Can I walk with the foot and heel?

No, you should not walk with this device on.

Can both products be used together?

Yes, the Devon foot and heel and the utility pad can be used at the same time to provide additional protection.

Can they be washed or cleaned?

No, washing or cleaning is not recommended. Dirty or damaged products should be disposed of and replaced.

Can they be cut?

Yes, the products can be cut to suit your clinical requirements, for positioning follow your Healthcare Professional's guidance.

This is a product overview and does NOT replace advice or guidance provided by your Healthcare Professional and the local NHS protocols.

If you have any further questions, please contact your Healthcare Professional on:

Name: _____

Tel: _____